

## PCR Protocol for Genotyping: PKD2<sup>neo</sup> (Pkd2cko)

### A. Digestion of mouse tail or ear notch, and **embryo tail (in red)**:

1. Add 100µL of Tissue Digestion Buffer and 2µL of Proteinase K per tail (~1-2mm length). For embryos tail add **50µL of Tissue Digestion Buffer and 1µL of Proteinase K**. Make sure tail is immersed in the buffer.
2. In a thermocycler incubate at 55°C for 1 h followed by 95°C for 8 min to inactivate the enzyme and hold at 10°C. For embryos incubate at 55°C for 30 min followed by 95°C for 8 min and hold at 10°C.
3. Vortex and store at 4°C (-20°C for long storage) or use immediately to set up the PCR.

### B. PCR Genotyping Protocol: PKD1<sup>neo</sup>

Primers			
<b>MG/5'flox11-13(A)-F</b>	5'-	<b>CCT TTC CTC TGT GTT CTG GGG AG</b>	-3'
<b>MG/5'flox11-13(B)-R</b>	5'-	<b>GTT TGA TGC TTA GCA GAT GAT GGC</b>	-3'

PCR Reaction		PCR Conditions			
BioMix (Bioline)	10.0 µL		Heated Lid		105°C
Primer F1	0.8 µL		Initial Denaturation	94°C	2 min
Primer B1	0.8 µL		Number of Cycles	x35	
Primer 2B	0.8 µL			94°C	30 sec
ddH <sub>2</sub> O	5.6 µL			56°C	30 sec
DNA template	2.0 µL			70°C	35 sec
Total Volume	20.0 µL		Final Extension	72°C	10 min
			Final Hold	10°C	

PCR Product Size (bp)	
Wild type band	232 bp
Pkd2 <sup>neo</sup> band	318 bp

### C. Reagents

Reagent	Cat #	Final Concentration	Working Concentration
<b>Tissue Digestion Buffer for ear notch or tail</b>			
Tris pH8.5		50mM	
EDTA		1mM	
Tween20		0.5%	
<b>Proteinase K (Invitrogen)</b>	25530-015	20mg/mL	
<b>BioMix (Bioline)</b>	BIO-25012		
<b>Primers</b>			10mM